

E-max™




Rapid Cadmium Testing in Food

Quantification of Unprecedented Low-level Cd By Monochromatic Excitation XRF

E-max is a monochromatic energy-dispersive x-ray fluorescence analyzer using a monochromatic beam optimized for Cd excitation. It delivers unprecedented Cd LOD taking advantage of the superior signal/background (S/B) ratio. It can rapidly survey and quantify low-level Cd in food with no or minimal sample preparation. Figure 1 shows Cd K-α signal peaks from three different rice samples.

Figure 1. Cd Signals for Different Rice Samples



Grain Application

A custom calibration curve for rice, as shown in Figure 3, is obtained using a set of reference rice standards. Various rice, wheat, and corn standard samples are validated based on this curve. The validation results are shown in Table 1 and Figure 4.

FP Approach Calibration

E-max applies fundamental-parameters (FP) approach to calibrate the system, using reference material to refine FP parameters. The FP approach allows E-max to measure various types of food with one calibration curve. Custom calibration curves can also be established to further improve accuracy for some specific sample types. Figure 2 shows Cd in dried rice, tomato and spinach.

Figure 2. Cd in Rice and Dried Vegetables




Figure 3. Cd in Rice





Table 1. Cd in Grain (Unit: ppm)

Sample	GBW10046 Wheat	GBW100496 Wheat	GBW08503C Wheat	GBW100493 Wheat	GBW100361 Rice	GBW100348 Rice	R080684a Rice	GBW100356 Rice	GBW100498 Corn	GBW100497 Corn
1	0.021	0.086	0.214	0.635	0.103	0.255	0.478	1.682	0.137	0.307
2	ND	0.105	0.204	0.648	0.097	0.237	0.452	1.741	0.108	0.296
3	0.027	0.098	0.223	0.638	0.132	0.232	0.482	1.691	0.102	0.280
Average	0.024	0.096	0.214	0.640	0.111	0.241	0.471	1.705	0.116	0.294
RSV	0.018	0.100	0.211	0.630	0.110	0.240	0.482	1.720	0.110	0.280
Accuracy	33.4%	3.6%	1.3%	1.6%	0.7%	0.6%	2.4%	0.9%	5.3%	5.1%

Figure 4 shows the validation results of rice, wheat, and corn standard samples .

Figure 4. Cd in Grain



Real Sample Testing

Cd in selected food samples in US markets has been tested with minimum sample preparation, as shown in Table 2.

Table 2. Cd in selected food samples (Unit: ppm)

Sample	Cereal	Chocolate	Dark Chocolate	Orzo Pasta	Salted Seaweed
					
1	ND	0.085	0.137	0.098	0.319
2	ND	0.077	0.131	0.099	0.328
3	0.019	0.089	0.143	0.109	0.331
Average	0.019	0.084	0.137	0.102	0.326

Repeatability and Accuracy

With the fundamental-parameter approach, E-max can be calibrated for various types of grain with high accuracy and excellent repeatability, as indicated in Table 3. The LOD of Cd is estimated to be 30ppb for 600s testing time.

Table 3. Repeated Measurements of Reference Standards and Others (Unit: ppm)

Sample	NMIJ7501	NMIJ7502	NIST1573	NIST1570	Pasta	Dark Chocolate
Testing time (s)	900	300	100	100	600	600
1	0.069	0.572	1.430	3.013	0.099	0.137
2	0.056	0.593	1.605	2.959	0.096	0.131
3	0.054	0.551	1.462	2.874	0.099	0.143
4	0.056	0.570	1.473	2.911	0.107	0.140
5	0.062	0.559	1.607	3.044	0.134	0.135
6	0.050	0.546	1.553	2.966	0.109	0.146
7	0.055	0.561	1.460	3.093	0.107	0.153
8	0.054	0.541	1.472	3.031	0.106	0.159
9	0.063	0.560	1.553	2.984	0.105	0.155
10	0.074	0.570	1.534	3.027	0.120	0.152
Average	0.059	0.562	1.515	2.990	0.108	0.145
RSV	0.052	0.548	1.517	2.876	NA	NA
SD	0.007	0.014	0.060	0.062	0.011	0.009
RSD	12.2%	2.5%	4.0%	2.1%	9.8%	6.1%
Accuracy	12.7%	2.5%	0.1%	3.8%	NA	NA

RSV: Reference Standard Value

SD: Standard Deviation

RSD: Relative Standard Deviation