

## **Oilseed Extractions (FOSFA comparable)**

XT10/XT10I and XT15/XT15I

Analytical Procedure Revised: 1/15/2020

The technique we use in the extraction of oilseeds is as follows:

- 1. Place 1.5 to 2.0 g of sample into the XT4 Filter Bag.
- 2. Dry the samples for three hours at 100° 105° C.
- 3. Remove the samples from the oven and place them in a desiccant pouch, allow them to cool and record their weights.
- 4. Extract the samples in an ANKOM XT instrument for forty minutes.
- 5. Dry the samples for thirty minutes at 100° 105° C.
- 6. Remove the samples from the oven and place them in a desiccant pouch, allow them to cool, and record their weights.
- 7. Place one XT4 Filter Bag (with sample) into a mortar and strike it with a pestle approximately 15-20 times being careful not to twist the pestle. There does not need to be a great deal of force when the pestle strikes the Filter Bag. Repeat for all Filter Bags.
- 8. Repeat steps 4-7 two more times (for a total of three times).
- 9. Add the three fat values together to get a combined fat total.